

Easy-to-use
BACTASLYDE®
Microbe Detection Device



Dip Slides



Pouch Products

RAKIRO BIOTECH SYSTEMS PVT LTD
WWW.BACTASLYDE.COM



BACTASLYDE : Comparison with Other Methods

Methods	No. of Steps involved	Trained manpower	Pre treatment of sample	Maintenance	Convenience	Capital Investment
Plate Count Method	6	Required	Not Required	Required	Tedious	Required
Most Probable Number	6	Required	Not Required	Required	Tedious	Required
Membrane Filter	6	Required	Not Required	Required	Tedious	Required
Direct Counts	4	Required	Required	Required	Convenient	Required
Bioluminescence	3	Required	Essential	Required	Convenient	Required
Direct Epifluorescence Technique (DEFT)	6	Required	Essential	Required	Convenient	Required
BACTASLYDE	1	Not Required	Not Required	Not Required	Most Convenient	Not Required

Easy-to-use

BACTASLYDE[®]

Microbe Detection Device

SALMONELLA BS PP2

Salmonella bacteria are zoonotic and can be transferred between humans and other animals. These are intestinal pathogens and cause many diseases viz. typhoid, paratyphoid, food poisoning etc. Classic typhoid fever is a serious and life-threatening infection which has a mortality approaching 20%. Cases of enteric fever almost always occur from the ingestion of contaminated food or water. These bacteria can survive for some time without a host for weeks outside a living body and they are not destroyed by freezing. They are frequently found in polluted water, contamination from the excrement of carrier animals being particularly important. This organism enters through the digestive tract due to ingestion of contaminated food. The "food infection" is generally thru' eating or handling contaminated raw or not properly cooked meat of cattle, swine, poultry, fish etc. FDA has laid down mandatory testing norms to control risk of salmonella contamination in foodstuff like spices, fish and eggs etc. BS PP2 is supplied in packs of 5 along with complete instructions for a qualitative analysis of the contamination.



How to use

Wash hands thoroughly with soap & water before handling BACTASLYDE BS-PP2



1. Media pouches & sterile bottle



2. Cut open both the pouches carefully



3. Pour the contents of both pouches into the bottle



4. Pour water upto the red arrow mark.



Observations

- Initially the sample alongwith the media will be a light green, grey or buff coloured turbid liquid.
- If Salmonella species are present, the liquid will turn dark green, dark brown or black.**

6. Then observe the colour of the liquid in the bottle



BACTASLYDE : Application in Different Industries

Industry	Application	Bactaslyde Code
Paint & Pigments	Process & production waters, raw materials, Finished products (water based paint & pigments)	BS 101, BS 102, BS 103
Metal Working Fluids	Cutting oils, coolants in the sump, raw water (make-up)	BS 101, BS 103, BS 115
Cooling Water Systems	Make-up waters, cooling waters, (recirculating water, basin water) Close systems	BS 103, BS 115, BS PP1 BS 125, BS 130
Food Processing Industry	Raw materials (e.g. cream, milk, meat, fish, vegetables), Water used for production and cleaning purpose, finished products (e.g. meat / fish preparations, cakes, soups, sauces, jams, squashes,spices)	BS 101, BS 102, BS 103, BS PP2, BS PP3, BS PP4
Pulp & Paper Industry	White waters (slurry and pulp mixtures), process waters	BS 101, BS 103, BS 115 BS 125
Sugar Industry	Primary, mixed, clarified juice, process waters, finished sugar	BS 101, BS 103, BS 115
Water Treatment Chemical Manufacturers	Laboratory and onsite, evaluation of biocides	BS 103, BS 115, BS PP 1 BS 125, BS 130

Industry	Application	Bactaslyde Code
Manufacturers Bulk Drug	Process waters, production water, (demineralized water)	BS 103, BS 115
Manufacturers Water & Waste Water Treatment Systems	Monitoring of waste water at different stages	BS 102, BS 103, BS 115 BS 125
Cosmetics	Process waters. (demineralized water) raw materials, finished product.	BS 101, BS 102, BS 103, BS PP 2, BS PP 3
Dairy	Raw milk, pasteurised milk, evaluation of the cleaning-in-process.	BS 101, BS 102, BS 103 , BS PP 2, BS PP 3
Brewery	Process waters, pasteurised Beer fermentation broth.	BS 101 , BS 102, BS 103, BS 115
Water Based Adhesives	Process waters, raw materials, finished product	BS 101, BS 103, BS 115
Oil and Petroleum	Injection waters, fuel (petrol, aviation)	BS 101, BS 103, BS 115
Fisheries	Ponds, Sea water, Processed products	BS 102, BS PP 2, BS PP 3, BS PP4

BS 101 - Yeasts & Fungi + TBC **BS 102** - Escherichia coli + TBC **BS 103** - Pseudomonas + TBC **BS 115** - SRB **BS 125** - Algae Species
BS PP1 - Iron Bacteria **BS PP 2** - Salmonella Species **BS PP 3** - Staphylococcus Species **BS PP 4** - Vibrio Species. **BS130** : Nitrifying / Denitrifying Bacteria

Our Other Product Range

AQUASOL

Analyzing Waters
Anytime... Anywhere...



Our range encompasses all known water parameters.

Acidity Alkalinity Amine Ammonium Arsenic Bromine
Calcium Hardness Carbon Dioxide Chloride Chlorine Chlorine Dioxide
Chloroscope Chromate Copper Dissolved Oxygen Fluoride
Hydrazine Hydrogen Peroxide Iron Molybdate Nitrate Nitrite
Orthophosphate pH Phosphonate Silica Sulphate Sulphite
Tannin Total Hardness Turbidity Zinc

ChekNsee

Power to see the Invisible Danger



Drinking water test kit

AQUASOL DIGITAL

*Electrochemistry Solutions for
all your applications with
a wide range of Digital Instruments.*



RAKIRO BIOTECH SYSTEMS PVT LTD

Navi Mumbai - 400 701, INDIA | Website: www.rakiro.net



An
ISO 9001:2015
Company