

**RAKIRO BIOTECH SYSTEMS PVT LTD**

R-466, TTC Industrial Area, MIDC
Rabale, Navi Mumbai - 400701
Tel No. +91-022-47804040
Email :- sales@rakiro.net

Doc No : TDSWT3010
Date : 01-02-2024
Type : UDAAQ
Product Code : WT3010

PRODUCT DATA SHEET**1 INFORMATION**

UDAAQ 3010-Instant Boiler Water Oxygen Scavenger

2 DESCRIPTION

UDAAQ 3010 is an instant oxygen scavenger. Product is recommended for the rapid and complete removal of dissolved oxygen from feed water and boiler water. The normal reaction time with oxygen is only 15 to 20 second at ambient temperatures.

3 SALIENT FEATURES

Instant Boiler Water Oxygen Scavenger
Ready to use liquid
Economical

Rapid action
Easy to handle and feed

4 ANALYTICAL DATA

Appearance	:	Pale yellow to reddish yellow color liquid
pH	:	Acidic
Density	:	1.0 To 1.3

5 HANDLING INSTRUCTIONS

Keep the container away from direct heat & sunlight. Keep the container closed when not in use. The product should not be swallowed and prolonged contact with the skin should be avoided. Should it come in contact with the eyes, flush with clean, cold water and get medical attention.

6 RECOMMENDED DOSAGE

An dosage of 25 PPM / DO is recommended but this dosage is to be adjusted in such a way that it should give protection to the system, while being economical. Dosage might vary depending on specific system parameter.