

**RAKIRO BIOTECH SYSTEMS PVT LTD**

R-466, TTC Industrial Area, MIDC

Rabale, Navi Mumbai - 400701

Tel No. +91-022-47804040

Email :- sales@rakiro.net

Doc No : TDSWT3030**Date :** 01-04-2024**Type :** UDAAQ**Product Code :** WT3030**PRODUCT DATA SHEET**> **INFORMATION**

Product Code: WT3030

Product Name: Oxygen Scavenger Hydrazine Based

Application: Boiler Treatment

> **DESCRIPTION**

The liquid oxygen scavenger, UDAAQ 3030 will instantly react with dissolved oxygen. Neither UDAAQ 3030 nor its reaction product will add to the solids content of the water. It also reduces ferric irons to the ferrous state, thus preventing harmful iron deposits within the boiler.

> **SALIENT FEATURES**

High Pressure Boiler Water Treatment

Environment-friendly and Non-Toxic Boiler Treatment

Economical and Available in ready to use Pack

High Temperatures and High-Pressure Stability

Easy to handle with Rapid Action for Instant solution

Complete Boiler Protection to enhance system life

> **ANALYTICAL DATA**

Appearance : COLOURLESS TO PALE YELLOW COLOUR LIQUID

pH : 9.00 - 13.00

Density : 0.95 - 1.05

> **HANDLING INSTRUCTIONS**

The use of protective Clothing, facemask, safety goggles, shoes and gloves are mandatory. Keep the container away from direct heat & sunlight. Keep the container closed when not in use. The product should not be swallowed and prolonged contact with the skin should be avoided. Should it come in contact with the eyes, flush with clean, cold water and get medical attention.

> **RECOMMENDED DOSAGE**

An dosage of 10 PPM / DO is recommended but this dosage is to be adjusted in such a way that it should give protection to the system, while being economical. Dosage might vary depending on specific system parameter.