

**RAKIRO BIOTECH SYSTEMS PVT LTD**

R-466, TTC Industrial Area, MIDC
Rabale, Navi Mumbai - 400701
Tel No. +91-022-47804040
Email :- sales@rakiro.net

Doc No : TDSWT3030
Date : 01-04-2024
Type : UDAAQ
Product Code : WT3030

PRODUCT DATA SHEET**1 INFORMATION**

Product Code: WT3030

Product Name: Oxygen Scavenger Hydrazine Based

Application: Boiler Treatment

2 DESCRIPTION

The liquid oxygen scavenger, UDAAQ 3030 will instantly react with dissolved oxygen. Neither UDAAQ 3030 nor its reaction product will add to the solids content of the water.

It also reduces ferric irons to the ferrous state, thus preventing harmful iron deposits within the boiler.

3 SALIENT FEATURES

Instant Boiler Water Oxygen Scavenger

Improves heat transfer

Imparts no solids to boiler water
solids to steam

Prevents iron oxide deposits

Stable

Imparts no Ready to use liquid

4 ANALYTICAL DATA

Appearance : Colourless to pale yellow colour liquid

pH : Alkaline

Density : 0.95 To 1.15

5 HANDLING INSTRUCTIONS

Keep the container away from direct heat & sunlight. Keep the container closed when not in use. The product should not be swallowed and prolonged contact with the skin should be avoided. Should it come in contact with the eyes, flush with clean, cold water and get medical attention.

6 RECOMMENDED DOSAGE

An dosage of 10 PPM / DO is recommended but this dosage is to be adjusted in such a way that it should give protection to the system, while being economical. Dosage might vary depending on specific system parameter.