

**RAKIRO BIOTECH SYSTEMS PVT LTD**

R-466, TTC Industrial Area, MIDC

Rabale, Navi Mumbai - 400701

Tel No. +91-022-47804040

Email :- sales@rakiro.net

**Doc No :** TDSWT8200**Date :** 01-04-2024**Type :** UDAAQ**Product Code :** WT8200**PRODUCT DATA SHEET**> **INFORMATION**

Product Code: WT8200

Product Name: Mill Sanitation Compound

Application: Sugar Processing

> **DESCRIPTION**

UDAAQ – 8200 is a chemical formulation consisting of synergistic blend of biocides which will sanitise the mill section area. UDAAQ 8200 restricts the inversion process of sugar in sugar juice. UDAAQ 8200 is a highly effective Microbiocide formulation wherein its chelating agent combines with the metal ions in the enzymes essential to the metabolism of microorganisms, thereby stunting their growth.

> **SALIENT FEATURES**

Environment-friendly and Non-Toxic Products

Economical and Ready to use in 5/20 KG Pack

Helps in Enhancing Sugar Extraction

High Chemical Stability with Enhanced Inhibition

Easy to handle with Rapid Action for Instant solution

Complete system Protection to Enhance life

> **ANALYTICAL DATA**

Appearance : PALE YELLOW TO AMBER COLOUR LIQUID

pH : 9.00 - 12.00

Density : 1.00 - 1.20

> **HANDLING INSTRUCTIONS**

The use of protective Clothing, facemask, safety goggles, shoes and gloves are mandatory. Keep the container away from direct heat & sunlight. Keep the container closed when not in use. The product should not be swallowed and prolonged contact with the skin should be avoided. Should it come in contact with the eyes, flush with clean, cold water and get medical attention.

> **RECOMMENDED DOSAGE**

Generally 5 – 10 ppm of UDAAQ 8200 is sufficient for mill sanitation. The dosage may vary w.r.t. stoppage of mills, cane quality, fiber percentage and quality of maceration water