

COD Digester

SKU: APCODDG1 | Range: 5 - 150 Degc

For easy and simple on site testing of COD (Chemical Oxygen Demand) and we have a specially designed COD Digester which helps in COD digestion process. Our COD reactor is used to digest closed micro reflux COD vial according to Standard Method and helps meet USEPA standards. It makes the process of COD Testing easier and safer than traditional open reflux titrimetric method.



KEY PRODUCT FEATURES

- Microprocessor based for fast and accurate measurement.
- Suitable for Standard Vial Size
- Additional Heating Programs available
- Predefined COD Program with Timer
- Auto Off and Alarm Functions
- Aluminum heating block with 12 vial capacity

TECHNICAL SPECIFICATIONS

| | |
|-------------------------------------|--|
| Parameter | |
| COD Digester / Dry Bath Incubator | |
| Model | |
| APCODDG1 | |
| Range | |
| 5 TO 150 DEG C | |
| Number of Blocks | |
| 12 Holes (Suitable for 16 MM Vials) | |
| Stability | |
| 0.5 DEG C | |
| Time Range | |
| 1 Min - 99H 59Min | |
| Power | |
| 240V AC 50 HZ (165W) | |
| Dimension | |
| 152 x 86 x 190 MM | |
| Meter Weight | |
| 2500 g | |